





FSCB Antibody

Product Code	CSB-PA708110LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5H9T9
Immunogen	Recombinant Human Fibrous sheath CABYR-binding protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibrous sheath CABYR-binding protein, FSCB, C14orf155
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FSCB


Image



Western blot
 All lanes: FSCB antibody at 1.5µg/ml + Mouse liver tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 88 kDa
 Observed band size: 88 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA708110LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA708110LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.