



Fxn Antibody, Biotin conjugated

Product Code	CSB-PA009086LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O35943
Immunogen	Recombinant Mouse Frataxin, mitochondrial protein (78-207AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Frataxin, mitochondrial (Fxn) (EC 1.16.3.1) [Cleaved into: Frataxin intermediate form; Frataxin mature form], Fxn, Frda
Immunogen Species	Mus musculus (Mouse)
Research Area	Tags & Cell Markers
Target Names	Fxn
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.