



FYCO1 Antibody, Biotin conjugated

Product Code	CSB-PA866262LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BQS8
Immunogen	Recombinant Human FYVE and coiled-coil domain-containing protein 1 protein (2-255AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	FYVE and coiled-coil domain-containing protein 1 (Zinc finger FYVE domain-containing protein 7), FYCO1, ZFYVE7
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	FYCO1