




FYTTD1 Antibody

Product Code	CSB-PA842744LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96QD9
Immunogen	Recombinant Human UAP56-interacting factor protein (1-318AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UAP56-interacting factor (Forty-two-three domain-containing protein 1) (Protein 40-2-3), FYTTD1, UIF
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FYTTD1

Image



Western blot

All lanes: FYTTD1 antibody at 4µg/ml

Lane 1: A375 whole cell lysate

Lane 2: Jurkat whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 36, 34, 29, 15 kDa

Observed band size: 36 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.