



FYTTD1 Antibody, Biotin conjugated

Product Code	CSB-PA842744LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96QD9
Immunogen	Recombinant Human UAP56-interacting factor protein (1-318AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UAP56-interacting factor (Forty-two-three domain-containing protein 1) (Protein 40-2-3), FYTTD1, UIF
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FYTTD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.