





GALK2 Antibody

Product Code	CSB-PA009199LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01415
Immunogen	Recombinant Human N-acetylgalactosamine kinase protein (161-458AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	N-acetylgalactosamine kinase (EC 2.7.1.157) (GalNAc kinase) (Galactokinase 2), GALK2, GK2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GALK2

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA009199LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreatic tissue using CSB-PA009199LA01HU at dilution of 1:100