



GBP5 Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA847203LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q96PP8 |
| Immunogen | Recombinant Human Guanylate-binding protein 5 protein (321-500AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Guanylate-binding protein 5 (EC 3.6.5.-) (GBP-TA antigen) (GTP-binding protein 5) (GBP-5) (Guanine nucleotide-binding protein 5), GBP5 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | GBP5 |