



GCSAM Antibody, Biotin conjugated

Product Code	CSB-PA854052LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N6F7
Immunogen	Recombinant Human Germinal center-associated signaling and motility protein (1-178AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Germinal center-associated signaling and motility protein (Germinal center B-cell-expressed transcript 2 protein) (Germinal center-associated lymphoma protein) (hGAL), GCSAM, GAL GCET2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	GCSAM
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.