



GHITM Antibody, Biotin conjugated

Product Code	CSB-PA872487LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H3K2
Immunogen	Recombinant Human Growth hormone-inducible transmembrane protein (293-345AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Growth hormone-inducible transmembrane protein (Dermal papilla-derived protein 2) (Mitochondrial morphology and cristae structure 1) (MICS1) (Transmembrane BAX inhibitor motif-containing protein 5), GHITM, DERP2 MICS1 TMBIM5
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	GHITM