



GK5 Antibody, FITC conjugated

Product Code	CSB-PA744414LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZS86
Immunogen	Recombinant Human Putative glycerol kinase 5 protein (96-301AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative glycerol kinase 5 (GK 5) (Glycerokinase 5) (EC 2.7.1.30) (ATP:glycerol 3-phosphotransferase 5), GK5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GK5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.