



GMEB1 Antibody, HRP conjugated

Product Code	CSB-PA896751LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y692
Immunogen	Recombinant Human Glucocorticoid modulatory element-binding protein 1 protein (2-245AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glucocorticoid modulatory element-binding protein 1 (GMEB-1) (DNA-binding protein p96PIF) (Parvovirus initiation factor p96) (PIF p96), GMEB1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GMEB1