



GPR146 Antibody, Biotin conjugated

Product Code	CSB-PA850271LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96CH1
Immunogen	Recombinant Human Probable G-protein coupled receptor 146 protein (301-333AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable G-protein coupled receptor 146 (G-protein coupled receptor PGR8), GPR146, PGR8
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GPR146
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.