



GPR17 Antibody, Biotin conjugated

Product Code	CSB-PA619758LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13304
Immunogen	Recombinant Human Uracil nucleotide/cysteinyl leukotriene receptor protein (1-36AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uracil nucleotide/cysteinyl leukotriene receptor (UDP/CysLT receptor) (G-protein coupled receptor 17) (P2Y-like receptor) (R12), GPR17
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GPR17
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.