




GSDMD Antibody

Product Code	CSB-PA009956LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P57764
Immunogen	Recombinant Human Gasdermin-D protein (1-256AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gasdermin-D (Gasdermin domain-containing protein 1) [Cleaved into: Gasdermin-D, N-terminal (GSDMD-NT); Gasdermin-D, C-terminal (GSDMD-CT)], GSDMD, DFNA5L GSDMDC1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	GSDMD

Image



Western blot

All lanes: GSDMD antibody at 4µg/ml

Lane 1: Jurkat whole cell lysate

Lane 2: HeLa whole cell lysate

Lane 3: A431 whole cell lysate

Lane 4: U251 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 53 kDa

Observed band size: 53 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.