



GSDMD Antibody, HRP conjugated

Product Code	CSB-PA009956LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P57764
Immunogen	Recombinant Human Gasdermin-D protein (1-256AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gasdermin-D (Gasdermin domain-containing protein 1) [Cleaved into: Gasdermin-D, N-terminal (GSDMD-NT); Gasdermin-D, C-terminal (GSDMD-CT)], GSDMD, DFNA5L GSDMDC1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	GSDMD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.