



GTF3C6 Antibody, FITC conjugated

Product Code	CSB-PA853368LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969F1
Immunogen	Recombinant Human General transcription factor 3C polypeptide 6 protein (2-213AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	General transcription factor 3C polypeptide 6 (Transcription factor IIIC 35 kDa subunit) (TFIIIC 35 kDa subunit) (TFIIIC35) (Transcription factor IIIC subunit 6), GTF3C6, C6orf51
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GTF3C6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.