



GTPBP6 Antibody, HRP conjugated

Product Code	CSB-PA04439B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43824
Immunogen	Recombinant Human Putative GTP-binding protein 6 protein (146-295AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative GTP-binding protein 6 (Pseudoautosomal GTP-binding protein-like), GTPBP6, PGPL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GTPBP6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.