



Gzma Antibody, HRP conjugated

Product Code	CSB-PA010081LB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11032
Immunogen	Recombinant Mouse Granzyme A protein (29-260AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Granzyme A (EC 3.4.21.78) (Autocrine thymic lymphoma granzyme-like serine protease) (CTLA-3) (Fragmentin-1) (T cell-specific serine protease 1) (TSP-1), Gzma, CtlA-3 CtlA3 Mtsp-1
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Gzma
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.