




HBD Antibody

Product Code	CSB-PA010152LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02042
Immunogen	Recombinant Human Hemoglobin subunit delta protein (2-147AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hemoglobin subunit delta (Delta-globin) (Hemoglobin delta chain), HBD
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	HBD

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA010152LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.