



HBE1 Antibody, Biotin conjugated

Product Code	CSB-PA010153LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02100
Immunogen	Recombinant Human Hemoglobin subunit epsilon protein (2-147AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hemoglobin subunit epsilon (Epsilon-globin) (Hemoglobin epsilon chain), HBE1, HBE
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	HBE1