



HBQ1 Antibody

Product Code	CSB-PA010159LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09105
Immunogen	Recombinant Human Hemoglobin subunit theta-1 protein (1-142AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hemoglobin subunit theta-1 (Hemoglobin theta-1 chain) (Theta-1-globin), HBQ1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	HBQ1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.