



# HDHD3 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA883592LC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9BSH5
<b>Immunogen</b>	Recombinant Human Haloacid dehalogenase-like hydrolase domain-containing protein 3 protein (1-251AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Haloacid dehalogenase-like hydrolase domain-containing protein 3, HDHD3, C9orf158
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	HDHD3
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.