



HEATR9 Antibody, Biotin conjugated

Product Code	CSB-PA003369LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A2RTY3
Immunogen	Recombinant Human Protein HEATR9 protein (1-219AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein HEATR9 (HEAT repeat-containing protein 9), HEATR9, C17orf66
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	HEATR9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.