



# HLA3 Antibody

Product Code	CSB-PA895005LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9XRX5
Immunogen	Recombinant Human HERV-H LTR-associating protein 3 protein (1-121AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	HERV-H LTR-associating protein 3 (Human endogenous retrovirus-H long terminal repeat-associating protein 3), HHLA3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	HLA3