



MFSD14B Antibody, HRP conjugated

Product Code	CSB-PA713104LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5SR56
Immunogen	Recombinant Human Hippocampus abundant transcript-like protein 1 protein (450-506AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hippocampus abundant transcript-like protein 1 (Major facilitator superfamily domain-containing 14B), MFSD14B, HIATL1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MFSD14B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.