



HIF1AN Antibody, Biotin conjugated

Product Code	CSB-PA865160LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NWT6
Immunogen	Recombinant Human Hypoxia-inducible factor 1-alpha inhibitor protein (2-254AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hypoxia-inducible factor 1-alpha inhibitor (EC 1.14.11.30) (EC 1.14.11.n4) (Factor inhibiting HIF-1) (FIH-1) (Hypoxia-inducible factor asparagine hydroxylase), HIF1AN, FIH1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HIF1AN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.