





HLA-DMA Antibody

Product Code	CSB-PA338892LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28067
Immunogen	Recombinant Human HLA class II histocompatibility antigen, DM alpha chain protein (27-233AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	HLA class II histocompatibility antigen, DM alpha chain (MHC class II antigen DMA) (Really interesting new gene 6 protein), HLA-DMA, DMA RING6
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	HLA-DMA

Image



Immunohistochemistry of paraffin-embedded human appendix tissue using CSB-PA338892LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA338892LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.