



HLA-DMA Antibody, FITC conjugated

Product Code	CSB-PA338892LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28067
Immunogen	Recombinant Human HLA class II histocompatibility antigen, DM alpha chain protein (27-233AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	HLA class II histocompatibility antigen, DM alpha chain (MHC class II antigen DMA) (Really interesting new gene 6 protein), HLA-DMA, DMA RING6
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	HLA-DMA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.