



HLA-DRB4 Antibody, Biotin conjugated

Product Code	CSB-PA03909D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13762
Immunogen	Recombinant Human HLA class II histocompatibility antigen, DR beta 4 chain protein (30-227AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	HLA class II histocompatibility antigen, DR beta 4 chain (MHC class II antigen DRB4), HLA-DRB4
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	HLA-DRB4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.