



# HLA-DRB4 Antibody, FITC conjugated

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA03909C0Rb  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                  |
| <b>Uniprot No.</b>         | P13762   |
| <b>Immunogen</b>           | Recombinant Human HLA class II histocompatibility antigen, DR beta 4 chain protein (30-227AA)  |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Human  |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | FITC   |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4               |
| <b>Purification Method</b> | >95%, Protein G purified   |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Alias</b>               | HLA class II histocompatibility antigen, DR beta 4 chain (MHC class II antigen DRB4), HLA-DRB4 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)   |
| <b>Research Area</b>       | Immunology   |
| <b>Target Names</b>        | HLA-DRB4   |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                    |