



# HMGCLL1 Antibody

<b>Product Code</b>	CSB-PA851531LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q8TB92
<b>Immunogen</b>	Recombinant Human 3-hydroxymethyl-3-methylglutaryl-CoA lyase, cytoplasmic protein (2-340AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	3-hydroxymethyl-3-methylglutaryl-CoA lyase, cytoplasmic (EC 4.1.3.4) (3-hydroxymethyl-3-methylglutaryl-CoA lyase-like protein 1) (Endoplasmic reticulum 3-hydroxymethyl-3-methylglutaryl-CoA lyase) (er-cHL), HMGCLL1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	HMGCLL1