



HMGCLL1 Antibody, Biotin conjugated

Product Code	CSB-PA851531LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TB92
Immunogen	Recombinant Human 3-hydroxymethyl-3-methylglutaryl-CoA lyase, cytoplasmic protein (2-340AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	3-hydroxymethyl-3-methylglutaryl-CoA lyase, cytoplasmic (EC 4.1.3.4) (3-hydroxymethyl-3-methylglutaryl-CoA lyase-like protein 1) (Endoplasmic reticulum 3-hydroxymethyl-3-methylglutaryl-CoA lyase) (er-cHL), HMGCLL1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HMGCLL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.