




HNRNPL Antibody

Product Code	CSB-PA010612LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14866
Immunogen	Recombinant Human Heterogeneous nuclear ribonucleoprotein L protein (89-335AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heterogeneous nuclear ribonucleoprotein L (hnRNP L), HNRNPL, HNRPL
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	HNRNPL

Image



Western Blot

Positive WB detected in: HepG2 whole cell lysate, Mouse heart tissue, Mouse brain tissue, Mouse kidney tissue


All lanes: HNRNPL antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 65, 51 kDa

Observed band size: 65 kDa



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA010612LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.