




HOPX Antibody

Product Code	CSB-PA859411LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BPY8
Immunogen	Recombinant Human Homeodomain-only protein (1-73AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeodomain-only protein (Lung cancer-associated Y protein) (Not expressed in choriocarcinoma protein 1) (Odd homeobox protein 1), HOPX, HOD HOP LAGY NECC1 OB1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	HOPX

Image



Immunohistochemistry of paraffin-embedded human heart tissue using CSB-PA859411LA01HU at dilution of 1:100