



HOPX Antibody, Biotin conjugated

Product Code	CSB-PA859411LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BPY8
Immunogen	Recombinant Human Homeodomain-only protein (1-73AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeodomain-only protein (Lung cancer-associated Y protein) (Not expressed in choriocarcinoma protein 1) (Odd homeobox protein 1), HOPX, HOD HOP LAGY NECC1 OB1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	HOPX