





HPCAL1 Antibody

Product Code	CSB-PA02269A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P37235
Immunogen	Recombinant Human Hippocalcin-like protein 1 protein (2-193AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hippocalcin-like protein 1 (Calcium-binding protein BDR-1) (HLP2) (Visinin-like protein 3) (VILIP-3), HPCAL1, BDR1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	HPCAL1


Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA02269A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA02269A0Rb at dilution of 1:100



Immunofluorescent analysis of PC-3 cells using CSB-PA02269A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.