



# HPCAL1 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA02269B0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P37235
<b>Immunogen</b>	Recombinant Human Hippocalcin-like protein 1 protein (2-193AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Hippocalcin-like protein 1 (Calcium-binding protein BDR-1) (HLP2) (Visinin-like protein 3) (VILIP-3), HPCAL1, BDR1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Neuroscience
<b>Target Names</b>	HPCAL1