



HSDL2 Antibody, Biotin conjugated

Product Code	CSB-PA761576LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6YN16
Immunogen	Recombinant Human Hydroxysteroid dehydrogenase-like protein 2 protein (1-345AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hydroxysteroid dehydrogenase-like protein 2 (EC 1.-.-.) (Short chain dehydrogenase/reductase family 13C member 1), HSDL2, C9orf99 SDR13C1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	HSDL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.