



Hspa1b Antibody, Biotin conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA010822HD01MO |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P17879 |
| Immunogen | Recombinant Mouse Heat shock 70 kDa protein 1B protein (2-642AA) |
| Raised In | Rabbit |
| Species Reactivity | Mouse |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Heat shock 70 kDa protein 1B (Heat shock 70 kDa protein 1) (HSP70.1), Hspa1b, Hcp70.1 Hsp70-1 Hsp70a1 Hspa1 |
| Immunogen Species | Mus musculus (Mouse) |
| Research Area | Others |
| Target Names | Hspa1b |