



HSPB7 Antibody, FITC conjugated

Product Code	CSB-PA10979C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBY9
Immunogen	Recombinant Human Heat shock protein beta-7 protein (1-170AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heat shock protein beta-7 (HspB7) (Cardiovascular heat shock protein) (cvHsp), HSPB7, CVHSP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HSPB7