



HSPB7 Antibody, FITC conjugated

Product Code	CSB-PA10979C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBY9
Immunogen	Recombinant Human Heat shock protein beta-7 protein (1-170AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heat shock protein beta-7 (HspB7) (Cardiovascular heat shock protein) (cvHsp), HSPB7, CVHSP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HSPB7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.