



HYDIN Antibody, HRP conjugated

Product Code	CSB-PA670104LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q4G0P3
Immunogen	Recombinant Human Hydrocephalus-inducing protein homolog protein (101-400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hydrocephalus-inducing protein homolog, HYDIN, HYDIN1 KIAA1864
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HYDIN