



ICA1L Antibody, FITC conjugated

Product Code	CSB-PA847691LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NDH6
Immunogen	Recombinant Human Islet cell autoantigen 1-like protein (1-187AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Islet cell autoantigen 1-like protein (Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 14 protein) (Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 15 protein), ICA1L, ALS2CR14 ALS2CR15
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ICA1L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.