



IDI2 Antibody, HRP conjugated

Product Code	CSB-PA871603LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXS1
Immunogen	Recombinant Human Isopentenyl-diphosphate delta-isomerase 2 protein (1-227AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Isopentenyl-diphosphate delta-isomerase 2 (EC 5.3.3.2) (Isopentenyl pyrophosphate isomerase 2) (IPP isomerase 2) (IPPI2), IDI2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	IDI2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.