



# IER2 Antibody, HRP conjugated

Product Code	CSB-PA01729B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BTL4
Immunogen	Recombinant Human Immediate early response gene 2 protein (1-223AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Immediate early response gene 2 protein (Protein ETR101), IER2, ETR101 PIP92
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	IER2