



IFIT5 Antibody, Biotin conjugated

Product Code	CSB-PA011023LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13325
Immunogen	Recombinant Human Interferon-induced protein with tetratricopeptide repeats 5 protein (107-320AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon-induced protein with tetratricopeptide repeats 5 (IFIT-5) (Interferon-induced 58 kDa protein) (Retinoic acid- and interferon-inducible 58 kDa protein) (P58), IFIT5, ISG58 RI58
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	IFIT5