



IGSF11 Antibody, Biotin conjugated

Product Code	CSB-PA691547LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5DX21
Immunogen	Recombinant Human Immunoglobulin superfamily member 11 protein (23-241AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Immunoglobulin superfamily member 11 (IgSF11) (Brain and testis-specific immunoglobulin superfamily protein) (Bt-IGSF) (V-set and immunoglobulin domain-containing protein 3), IGSF11, BTIGSF CXADRL1 VSIG3
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IGSF11