



IGSF11 Antibody, HRP conjugated

Product Code	CSB-PA691547LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5DX21
Immunogen	Recombinant Human Immunoglobulin superfamily member 11 protein (23-241AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Immunoglobulin superfamily member 11 (IgSF11) (Brain and testis-specific immunoglobulin superfamily protein) (Bt-IGSF) (V-set and immunoglobulin domain-containing protein 3), IGSF11, BTIGSF CXADRL1 VSIG3
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IGSF11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.