



IMPDH2 Antibody, Biotin conjugated

Product Code	CSB-PA011701LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12268
Immunogen	Recombinant Human Inosine-5'-monophosphate dehydrogenase 2 protein (5-514AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inosine-5'-monophosphate dehydrogenase 2 (IMP dehydrogenase 2) (IMPD 2) (IMPDH 2) (EC 1.1.1.205) (IMPDH-II), IMPDH2, IMPD2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	IMPDH2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.