



IPO11 Antibody, Biotin conjugated

Product Code	CSB-PA887057LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UI26
Immunogen	Recombinant Human Importin-11 protein (1-270AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Importin-11 (Imp11) (Ran-binding protein 11) (RanBP11), IPO11, RANBP11
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	IPO11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.